


▶ Reflexive Evaluation of Serum Protein Electrophoresis

In order to maximise the clinical utility of the Serum Protein Electrophoresis (SPE), in a cost-effective manner, RDL routinely, and at no additional charge, perform both Immunofixation Electrophoresis (IFE) and Quantitative Immunoglobulins (QIG) on sera suspected of having a monoclonal spike. Furthermore, QIG will be determined if the SPE suggests hypogammaglobulinemia.

▶ Summary of Algorithm



RDL's SPEP, a reflexive, cost-effective evaluation of your patients with suspected or known inflammatory or neoplastic disorders. Specimen requirement is 2 mL serum, ambient, refrigerated or frozen.